

GSure® Fungal DNA kit

#G45331 50
preparations

Store at Room Temperature

Procedure:

1. Take ~25-50mg (wet weight) of fungi in a mortar, add 250 µl GDFU1 buffer to it. Add 250 µl either Part A or Part B. For most of the fungi, use part A. Part B is to be used when working with fungi species that produces high amount of secondary metabolites.
2. Collect 300µl-400µl of slurry in a fresh microfuge tube and vortex vigorously. Fungi particle should also come with the buffer.
3. Incubate the tube at 70°C for 30 minutes and vortex after every 5 minutes. Color of the suspended cells may turn dark green at this time.
4. After incubation step, centrifuge the tube at maximum speed (10,000xg) for 10 minutes. Collect 250 µl of the clear supernatant in a fresh microfuge tube,
5. Add 250 µl Buffer GDFU2 and mix by inverting the tube 4–6 times.
6. Add 350 µl Buffer GDFU3 and invert the tube immediately. Shake vigorously to mix the solutions, **DO NOT VORTEX AT THIS STAGE.** Vortex may cause shearing of genomic DNA.
7. Apply the mixture to the GMini Spin Column by decanting or pipetting. **Avoid mixing of cell debris with the supernatant** as this may clog GMini spin Column thus lowering the DNA yield.
8. Centrifuge at 8500xg for 30–60 s. Discard the flow-through.
9. Wash GMini Spin Column by adding 600µl Membrane Wash Buffer* and centrifuging for 30–60 s as previously.
10. Discard the flow-through.
11. Repeat washing step.

12. Discard the flow-through, and centrifuge for an additional 2 min to remove residual wash buffer from membrane.

This step is extremely important to ensure complete removal of ethanol. Presence of ethanol in purified DNA may inhibit successive enzymatic reactions.

13. Place the GMini Spin Column in a clean 1.5 ml microcentrifuge tube (not provided). To elute DNA, add 50 µl Nuclease-free water to the center of each GMini Spin Column, let stand for 1 min, and centrifuge for 1 min at maximum speed (~8500Xg) on a table top microcentrifuge.
14. Discard the column and collect the eluted DNA present in microcentrifuge tube.
 - If required, increase the volume of GDFU1, Part A or Part B, GDFU2 and GDFU3 accordingly.
 - **If any sediment found in any of the isolation buffers, warm the containers at 50°C until it dissolves.**

Kit Contents:

GDFU1 Buffer	: 15ml
Part A	: 15ml
Part B	: 15ml
GDFU2 Buffer	: 15ml
GDFU3 Buffer	: 20ml
*Membrane Wash Buffer	: 30ml
Spin Column	: 50 pcs.
Nuclease-free Water	: 3ml

Reconstitution of Membrane Wash Buffer:

Before using the kit for first time, add 45ml of absolute ethanol (molecular biology grade) with the provided wash buffer. Mix thoroughly by shaking. Once ethanol is added, tighten the cap properly after each use.