

## 2X Prime Taq PCR Master Mix

#G7118A                      250 Reactions

Store at -20°C

Spin tubes briefly before use

| Catalogue Number | Pack Size     | ✓ |
|------------------|---------------|---|
| G7118            | 100 Reaction  |   |
| G7118A           | 250 Reaction  | ✓ |
| G7118B           | 1000 Reaction |   |

### Introduction

Prime Taq DNA Polymerase is an optimized combination of Taq DNA polymerase and high fidelity DNA polymerases from Pyrococcus species for use in routine and difficult PCR experiments. The 3'→5' exonuclease activity of the high fidelity DNA Polymerase increases the fidelity and robustness in the amplification by Taq DNA Polymerase, even from very low copy number of template. Taq DNA polymerase is heat-stable and will synthesize DNA at elevated temperatures from single-stranded templates in the presence of primer. The formulation of this enzyme mix has been optimized such a way to improve its yield and activity for long PCR. It can be used for PCR of up to 10 kb product.

### Guidelines for using Hi-Prime Taq DNA Polymerase

The Master Mix contains Prime Taq DNA reaction buffer, 5mM MgCl<sub>2</sub>, and 0.4mM of each dNTP.

### Package information

|                              |   |
|------------------------------|---|
| <b>G7118 100 Reactions</b>   | Material provided: 2X PCR Master Mix (2 x 1250 ul), 10ul Control DNA template, 10ul Control primer Mix        |
| <b>G7118A 250 Reactions</b>  | <b>Material provided: 2X PCR Master Mix (5 x 1250 ul), 10ul Control DNA template, 10ul Control primer Mix</b> |
| <b>G7118B 1000 Reactions</b> | Material provided: 2X PCR Master Mix (20 x 1250 ul), 10ul Control DNA template, 10ul Control primer Mix       |

### Unit definition

One unit incorporates 10nmol of deoxy-ribonucleotide into acid-insoluble product in 30 minutes at 74°C. Unit assay conditions: 25 mM TAPS (pH 9.3), 50 mM KCl, 2 mM MgCl<sub>2</sub>, 1 mM DTT, 0.2 mM dATP, dCTP, dGTP, dTTP utilizing M13mp18DNA as template.

### Enzyme

The amount of enzyme will be 1.25 Unit in 50 µl reactions finally. When cloning fragments amplified with Prime Taq DNA Polymerase a extra adenine overhang base is incorporated at the 3' end of the amplified PCR product. That A overhang is required to perform the TA cloning.

**Table 1. Pipetting instructions** (add items in this order).

| Component         | 50 reaction | Final conc. |
|-------------------|-------------|-------------|
| H <sub>2</sub> O  | To 50µl     |             |
| 2x PCR Master Mix | 25 µl       | 1X          |
| dNTPs             | -           | 200µM each  |
| MgCl <sub>2</sub> | -           | 2.5mM       |
| Forward primer    | X µl        | 0.5µM       |
| Reverse Primer    | X µl        | 0.5µM       |
| Template DNA      | X µl        |             |

**Table 2. Cycling Instructions**

| Cycle step           | 2-step Protocol |           | Cycles |
|----------------------|-----------------|-----------|--------|
|                      | Temp.           | Time      |        |
| Initial Denaturation | 94°C            | 5min      | 1      |
| Denaturation         | 94°C            | 15-30sec  | 25-35  |
| Annealing            | 45-63°C         | 15-60sec  |        |
| Extension            | 72°C            | 60sec/1kb |        |
| Final extension      | 72°C            | 10 min    | 1      |
| Store                | 4°C             | Hold      |        |

### Template

For low complexity DNA (e.g. plasmid, lambda or BAC DNA): 2 pg– 20 ng per 50 µl reaction volume needed to be used. For high complexity genomic DNA, the amount of DNA template should be 20–200 ng per 50 µl reaction volume. If cDNA synthesis reaction mixture is used directly as a source for the template, the volume used should not exceed 10 % of the final PCR reaction volume.

### Primers

The recommendation for final primer concentration is 0.5 µM. If required, the primer concentration may be optimized between 0.2–1.0 µM. The results from primer T<sub>m</sub> calculations can vary significantly depending on the method used. Always use the T<sub>m</sub> calculator and instructions from reputed website to determine the T<sub>m</sub> values of primers and optimal annealing temperature. If using a two-step PCR protocol, where both primer annealing and extension occur in a single step at 72°C, the primers should be designed accordingly.

### Note

- Use 94°C for Denaturation
- The annealing temp depends on the primer.
- Use extension at 72°C for upto 6kb amplification and at 68°C for above 6kb amplification.
- Use 200 µM of each dNTP. Do not use dUTP.
- Note: PrimeTaq DNA Polymerase produces A overhang at the 3' end.

### Mg<sup>2+</sup> and dNTP

2X G9 Taq PCR Mater Mix is optimized to work well for most amplicons. Standard concentration of MgCl<sub>2</sub> in final reaction condition is 2.5 mM.

### DMSO

The recommended reaction conditions for GC-rich templates include 3 % DMSO as a PCR additive, which aids in the denaturing of templates with high GC contents. For further optimization the amount of DMSO should be increased in 2 % increments. In some cases DMSO may also be required for supercoiled plasmids to relax for denaturation. Other PCR additives such as formamide, glycerol, and betaine are also compatible with G9 Taq DNA Polymerase. If high DMSO concentration is used, the annealing temperature must be decreased, as DMSO affects the melting point of the primers. It has been reported that 10 % DMSO decreases annealing temperature by 5.5–6.0°C.

### Troubleshooting

| Optimize annealing temperature.   | Increase extension time.   |
|---|--|
| <ul style="list-style-type: none"> <li>• Optimize annealing temperature.</li> <li>• Repeat and make sure that there are no pipetting errors.</li> </ul> | <ul style="list-style-type: none"> <li>• Increase the number of cycles.</li> </ul>   |
| <ul style="list-style-type: none"> <li>• Make sure the cycling protocol was performed as recommended.</li> </ul>  | <ul style="list-style-type: none"> <li>• Titrate DMSO (2–8 %) in the reaction.</li> </ul>  |
| <ul style="list-style-type: none"> <li>• Use fresh high-quality dNTPs.</li> </ul>   | <ul style="list-style-type: none"> <li>• Denaturation temperature may be too low. Optimal denaturation temperature for most templates is 94–98°C.</li> </ul> |
| <ul style="list-style-type: none"> <li>• Do not use dNTP mix or primers that contain dUTP or dNTP.</li> </ul>   | <ul style="list-style-type: none"> <li>• Optimize denaturation time.</li> </ul>  |
| <ul style="list-style-type: none"> <li>• Titrate template amount. Template DNA may be damaged. Use freshly prepared template.</li> </ul>                | <ul style="list-style-type: none"> <li>• Check the purity and concentration of the primers.</li> <li>• Check primer design.</li> </ul>                       |
| Non-specific products - High molecular weight smears  |  |
| <ul style="list-style-type: none"> <li>• Make sure the extension time used was not too long.</li> </ul>   | <ul style="list-style-type: none"> <li>• Titrate template amount.</li> </ul>   |

- Reduce the total number of cycles.
- Increase annealing temperature or try 2-step PCR protocol.
- Optimize denaturation temperature.
- Decrease primer concentration.

### Non-specific products - Low molecular weight discrete bands

- Increase annealing temperature
- Make sure the extension time used was not too long.
- Titrate template amount.
- Decrease primer concentration.
- Design new primers.

| Component                | Volume |
|--------------------------|--------|
| 2X G9 Taq PCR Master Mix | 25 µl  |
| Control DNA template     | 1 µl   |
| Control primer Mix       | 1 µl   |
| Nuclease-free water      | 23 µl  |

### Control PCR cycling conditions

| Cycling conditions |               |
|--------------------|---------------|
| Step 1             | 94°C – 5 min  |
| Step 2             | 94°C – 30s    |
| Step 3             | 57°C – 30s    |
| Step 4             | 72°C – 4 min  |
| Step 5             | 72°C – 10 min |
| Step 6             | 4°C – Hold    |

Repeat Steps 2 to 4 for 25 Cycles. Control PCR will produce a 2Kb product.